

SurgiTel®

SurgiCam® PRO

Next Generation Digital Video Camera

Easy to Use

Compact

Versatile



NEW

"...it allows me to show
and record exactly what
I'm seeing!"

Dr. R. Tyler Green

A Breakthrough in Medical Technology


The SurgiCam® Pro Digital Video Camera is the industry's first digital surgical camera able to capture DVD-quality video directly via USB for digital recording on a desktop or laptop computer.

This patented technology offers a unique way to visually document clinical procedures in digital format, thus eliminating the need for additional specialized control units.

- USB plug & play technology
- Easy to use interface
- Advanced miniaturization
- Superior image sensor and lens
- High resolution snapshot option
- Audio recording option



Small and light weight



Easy to use interface

| SurgiCam® Pro Features | |
|---------------------------------|--------------------------------------------------------------------------------------------------------------------------|
| Length (with 3x lens)/ Diameter | 40mm (with audio) / 22mm 36mm (without audio) / 22mm |
| Weight | 13.5g (with audio) 10.5g (without audio) |
| USB Cord | Length: up to 14 feet EMI Filter: built in |
| System Requirements | - Windows XP, Vista or 7 - 1.2GHz or faster processor - 512MB minimum system RAM - 50GB minimum hard disk space |
| Lens Options | 2x, 3x and 5x |
| Audio Recording | Yes (optional) |
| Video Resolution | 640x480 pixels. 30 frames/sec. , AVI |
| Snapshot Resolution | 1280x1024, JPEG |
| Warranty | 2 years |



Optional computer stand